Vitrex Sterilance Press II

Sterile Safety Lancet

Technical Data Sheet



PRODUCT:


The Vitrex safety lancets are used to carry out skin punctures, in order to collect micro-samples of capillary blood.

- Sterile and disposable.
- Safety lancets with activation button.
- Ergonomic design case, created to provide a safe environment.
- Easy to use: easy use, 2-step operation and no charging.
- Safety, after use the lancet retracts inside the case, avoiding accidents.

The dimensions:

- Length: 39.3 ±1.5 mm
- Width: 14.8 ±0.5 mm

t



Needle core
Spring
Steel needle
Trigger
End cap
Shell
Protection cap

MATERIALS:

The materials used to manufacture our Vitrex safety lancets are Medical grade:

Component	Material
Needle	Stainless Steel (06Cr19Ni10)
Shell	ABS – QS-0613-017
Button	HDPE-QS-0613-017

PRINT:

The terminology, graphic symbols and information applied to this medical device comply with EN 1041:2013 och ISO 15223-1:2016.

QUALITY:

- Quality Management System in acccordance with ISO 9001:2015 and ISO 13485:2016.
- Complies with Annex II, (Excluding Section 4) of the Directive 93/42/EEC.
- Sterilization by gamma radiation.

CHARACTERISTICS

- Shelf life-time: 5 years
- Disposable

TECHNICAL SPECIFICATIONS:

Colour Code		Configuration	Penetration Depth	Article No.	Application
Green		18G blade	1.8 mm	393213	Special baby care and maternity units for heel puncture
Blue		21G needle	1.8 mm	397213	As Yellow but for increased blood flow
Orange		21G needle	2.2 mm	395213	Haemoglobin (blood banks), cholesterol, blood typing, coagulation, blood gas and other diagnostic tests
Pink		21G needle	2.8 mm	396213	As Orange but for thicker adult skin
Yellow		26G needle	1.8 mm	394213	Diabetes clinics/hospital wards for blood glucose monitoring
Purple		28G needle	1.8 mm	398213	Low blood flow, less pain
Turquoise		30G needle	1.8 mm	399213	Low blood flow, less pain

Shelf life 5 years

Vitrex Medical a/s
Vasekaer 6-8 - DK 2730 Herlev
Phone +45 4494 7011
vitrex@vitrex.dk
www.vitrexmedical.com

Vitrex Sterilance Press II

Sterile Safety Lancet

Technical Data Sheet



PACKAGING:

The inner box is made of white cardboard and the outer carton is made of corrugated cardboard.
Each box contains batch number, expiration date and sterilization information.

- An inner box contains 100 individual Vitrex safety lancet units.
- One outer carton contains 10 inner boxes.

Inner Box

Length	Width	Height	Quantity
13.5 cm	9 cm	7.5 cm	100 pcs



Transit packaging Pallet

Packaging:

Primary packaging:	100 pcs
Transit packaging:	1,000 pcs
Pallet:	98,000 pcs

LANGUAGES:

The labels contain 13 languages and comply accordingly to the European directive and may be sold in countries where one of these languages are accepted by national law.:

- Danish
- Spanish
- Italian
- Arabic
- Greek
- Swedish
- German
- Portuguese
- Finnish
- English
- French
- Russian
- Dutch

Vitrex Medical a/s
Vasekaer 6-8 - DK 2730 Herlev
Phone +45 4494 7011
vitrex@vitrex.dk
www.vitrexmedical.com