Product Specification Sheet

VWR MAGNETIC STIRRING BARS, CYLINDRICAL, WITH REMOVABLE PIVOT RING

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionés respectent toutes les specifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Length	Ø	Pk
442-0369	12 mm	8 mm	10
442-0371	25 mm	8 mm	10
442-0372	25 mm	10 mm	10
442-0373	32 mm	8 mm	10
442-0374	38 mm	8 mm	10
442-0375	38 mm	10 mm	10
442-0376	42 mm	10 mm	5
442-0377	45 mm	8 mm	5
442-0378	51 mm	8 mm	5
442-0379	51 mm	10 mm	5
442-0381	75 mm	12 mm	5
442-0382	102 mm	16 mm	1
442-0383	127 mm	16 mm	1
442-0384	150 mm	19 mm	1

Material: PFTE covered, strong Alnico V magnetic core

Properties: Good chemical resistance

Temperature Resistance: -200° up to +250°C

Sterilisation: Chemical and thermal, do not gamma irradiate

Compliant to: FDA

Product Specification Sheet

VWR MAGNETIC STIRRING BARS, CYLINDRICAL, WITH REMOVABLE PIVOT RING

PTFE used in the manufacture these products are pure, virgin material which conforms to FDA and USP class VI requirements and the following specified substance content. It does not contain, or use in its manufacture any Class 1 ODC according to the Clean Air Act 1990, PFOA, POP, DEHIP and BPA, no animal material is used in the manufacture or is present in this material.

	Description	Content (%)
Group A	Cadmium and Cadmium compounds	Not detected
	Hexavalent Chromium compounds	Not detected
	Lead and Lead compounds	Not detected
	Mercury and Mercury compounds	Not detected
	Tributyl Tin Oxide (TBTO)	Not detected
	Tributyl Tins & Triphenyl Tins	Not detected
Group B	Polybrominated Biphenyls (PBBs)	Not detected
	Polybrominated Diphenyl ethers (PBDE)	Not detected
	Polychlorinated Biphenyls (PCBs)	Not detected
	Polychloronaphthalenes	Not detected
	Short Chain Chlorinated Paraffins	Not detected
Group C	Asbestos	Not detected
	Azo Colorants	Not detected
	Ozone Depleting Substances	Not detected
	Radioactive Substances	Not detected

Disclaimer: *VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.*

Date of issue: **October 2012**