Product Specification Sheet

VWR SPECIMEN CONTAINERS, PP

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release.The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionés respectent toutes les specifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No	Volume (ml)	Shape	Cap	Sterile	ø lower (mm)	ø higher (mm)	Height (mm)	Unit/Inner Package	Pack Size
216-0855	120	Conical	Separate red screw cap	N	47,9	56	71,5	50	500
216-0856	120	Conical	Red screw cap	Y	47,9	56	71,8	Individually	250
216-0857	120	Conical	Blue screw cap	Y	47,9	56	71,8	/	300
216-0858	120	Conical	Blue screw cap	N	47,9	56	71,8	/	300
216-0859	120	Conical	Separate red screw cap	N	47,9	56	71,5	/	500
216-0860	120	Conical	White screw cap	Y	47,9	56	71,8	/	300
216-0861	120	Conical	Red screw cap	Y	47,9	56	71,8	/	300
216-0862	120	Conical	Green Screw cap	Y	47,9	56	71,8	/	300
216-0863	120	Conical	White screw cap	N	47,9	56	71,8	/	300
216-0864	120	Conical	Red screw cap	N	47,9	56	71,8	/	300
216-0865	120	Conical	Green screw cap	N	47,9	56	71,8	/	300
216-0866	60	Cylindrical	Red screw cap & spatula	N	37,3	38,2	64,3	/	600
216-0867	60	Cylindrical	Red screw cap & spatula	Y	37,3	38,2	64,3	/	600
216-0868	120	Conical	Red screw cap	N	47,9	56	71,8	Individually	250
216-0869	120	Conical	White screw cap	Y	47,9	56	71,8	Individually	250

Intended use: Collection of human specimens, other than blood.
For single use

Container Material: Polypropylene

Cap Material: Polyethylene

VWR International bvba/sprl, Haasrode Researchpark 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium
http://www.vwr.com
Version: 2.0

Product Specification Sheet

VWR SPECIMEN CONTAINERS, PP

Compliant: with the essential requirements of IVD Directive 98/79/EC.
tested acc. To EN 14254

Properties: Transparent, resistant at temperature between 0° C to 120° C.
Do not expose or leave near open flames

Chemical Properties: Excellent chemical resistance to acids, alkali and oils.

Shelf life: 5 years

Disclaimer: *VWRwants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there's any doubt ask us.*

Date of issue: **June, 2012**

VWR International bvba/sprl, Haasrode Researchpark 2020, Geldenaaksebaan 464, 3001 Leuven, Belgium
http://www.vwr.com