BRAND GMBH + CO KG
Postfach 11 55
97861 Wertheim
Germany
Tel: +49 9342 808-0
Fax: +49 9342 808-98000
info@brand.de
www.brand.de

BRAND

BRAND GMBH + CO KG · Otto-Schott-Straße 25 · 97877 Wertheim · Germany

Wertheim, 05/09/2018

Product information

Beaker, low form, PP, blue printed scale
Cat. No.: 87608

The above mentioned article is made of PP, high clarity. ISO 7056. With graduation, beaded rim and spout.
Cleaning up to max. 60 °C is recommended to preserve marks and inscriptions.

Capacity ml	Subdivision ml	Height mm	max. Ø mm	Pack of
100	10	70	55	12

For further information please have a look at our General Catalog 900 page 266 or visit our online catalog.

With best regards,

BRAND GMBH + CO KG
i.A.

Lena Scheuermann
- Customer Service-