DURAN
WHEATON
KIMBLE

Datenblatt - *Data Sheet*

Article	**Test Tube 125x16x0,6-0,7** Boro 5.1, without rim round bottom
Article Number	**90105043** (Alternative: A512516000711)
Material:	**Borosilicate Glass 5.1** •Water resistance tested to ISO 719 (at 98 °C): HGB 1 •Alkali resistance tested to ISO 695: class A2 •Acid resistance tested to ISO 12116: S1 •Maximum using temperature long-time load 400°C (*) •Maximum using temperature Short-time load 500°C (*) •Autoclavable
Dimension/Standard:	**Gemäß ISO 4142 / refering to ISO 4142**

Over all length H	Outside diameter D	Wall W
125 mm	16 mm	0,7 mm

