Technical Data Sheet



VWR® 15mL High-Performance Centrifuge Tubes

Description: 15mL High-Performance Centrifuge Tubes with Flat or Plug Caps, Sterilized/Non-Sterilized

Purpose: Centrifuge tubes are popular in samples centrifuging or storage in many research area.


Materials


Tube: PP (Polypropylene), Compliance with USP Class VI Color: Clear

Cap: HDPE (High Density Polyethylene), Compliance with USP Class VI Color: Green

Rack: PP (Polypropylene), Compliance with USP Class VI Color: Green

Dimensions (Unit:mm)





15mL Centrifuge Tubes with Plastic Rack

Features

- Easy-to-read black graduations in ±2% increments
- Contains a large, white frosted writing area
- Autoclavable at 121℃ and freezable to -80℃
- Meet the requirements of BSE/TSE
- Latex Statement: This product is latex free
- Non-pyrogenic: The acceptance level for product is 0.05 EU/ml. (Tachypleus Amebocyte Lysate (TAL) Clot-gel test)
- Meet the requirements of RoHS
- Leak-proof
- DNase/RNase free
- Rack marked with numbers (1-25) for identification
- Angular at one end to stack up easily
- Rack can be washed and reused
- Sterilized by gamma irradiation SAL10^{-6} (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system

② STERILE R

VWR® 15mL High-Performance Centrifuge Tubes

VWR EU EAN. No.	VWR NA Cat. No.	Capacity (mL)	Maximum RCF	Tube Shape	Tube Color	Cap Shape	Cap Color	Sterilization	Packaging Configuration
525-1068	89039-664	15mL	12,500g	Conical-Bottom	Clear	Flat cap	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1070	89039-668	15mL	12,500g	Conical-Bottom	Clear	Plug cap	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1069	89039-666	15mL	12,500g	Conical-Bottom	Clear	Flat cap	Green	Yes	25/RK, 2 RK/PK, 10 PK/CS, 500/CS
525-1071	89039-670	15mL	12,500g	Conical-Bottom	Clear	Plug cap	Green	Yes	25/RK, 2 RK/PK, 10 PK/CS, 500/CS
525-1066	89401-566	15mL	12,500g	Conical-Bottom	Clear	Flat cap	Green	No	Tube:500/Bag, Cap:500/Bag, 500/CS
525-1067	89401-568	15mL	12,500g	Conical-Bottom	Clear	Plug cap	Green	No	Tube:500/Bag, Cap:500/Bag, 500/CS

Technical Data Sheet



VWR® 50mL High-Performance Centrifuge Tubes

Description: 50mL High-Performance Centrifuge Tubes with Flat or Plug Caps, Sterilized/Non-Sterilized

Purpose: Centrifuge tubes are popular in samples centrifuging or storage in many research area.


Materials

Tube: PP (Polypropylene), Compliance with USP Class VI Color: Clear


Cap: HDPE (High Density Polyethylene), Compliance with USP Class VI Color: Green

Rack: PP (Polypropylene), Compliance with USP Class VI Color: Green


Dimensions (Unit:mm)



50mL Centrifuge Tube Plug Cap Flat Cap



50mL Freestanding Centrifuge Tube Plug Cap



50mL Centrifuge Tubes with Plastic Rack

Features

- Easy-to-read black graduations in ±2% increments
- Contains a large, white frosted writing area
- Autoclavable at 121℃ and freezable to -80℃
- Meet the requirements of BSE/TSE
- Latex Statement: This product is latex free
- Non-pyrogenic: The acceptance level for product is 0.05 EU/ml. (Tachypleus Amebocyte Lysate (TAL) Clot-gel test)
- Meet the requirements of RoHS
- Leak-proof
- DNase/RNase free
- Rack marked with numbers (1-25) for identification
- Angular at one end to stack up easily
- Rack can be washed and reused
- Sterilized by gamma irradiation SAL10^{-6} (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® 50mL High-Performance Centrifuge Tubes

VWR EU EAN. No.	VWR NA Cat. No.	Capacity (mL)	Maximum RCF	Tube Shape	Tube Color	Cap Shape	Cap Color	Sterilization	Packaging Configuration
525-1075	89039-656	50mL	12,500g	Conical-Bottom	Clear	Flat cap	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1077	89039-660	50mL	12,500g	Conical-Bottom	Clear	Plug cap	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1074	89039-658	50mL	12,500g	Conical-Bottom	Clear	Flat cap	Green	Yes	25/RK/PK, 20 PK/CS, 500/CS
525-1076	89039-662	50mL	12,500g	Conical-Bottom	Clear	Plug cap	Green	Yes	25/RK/PK, 20 PK/CS, 500/CS
525-1072	89401-562	50mL	12,500g	Conical-Bottom	Clear	Flat cap	Green	No	Tube:500/Bag, Cap:500/Bag, 500/CS
525-1073	89401-564	50mL	12,500g	Conical-Bottom	Clear	Plug cap	Green	No	Tube:500/Bag, Cap:500/Bag, 500/CS
525-1078	89079-494	50mL	12,500g	Freestanding	Clear	Plug cap	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1079	89079-492	50mL	12,500g	Freestanding	Clear	Plug cap	Green	No	50/Bag, 10 Bags/CS, 500/CS
525-1080	89079-490	50mL	12,500g	Freestanding	Clear	Plug cap	Green	No	Tube:500/Bag, Cap:500/Bag, 500/CS

VWR®

Technical Data Sheet



VWR® 15mL Ultra High Performance Centrifuge Tubes


Description: 15mL Ultra High Performance Centrifuge Tubes with Flat or Plug Caps, Flat cap with sealing ring, Sterilized/Non-Sterilized

Purpose: Centrifuge tubes are popular in samples centrifuging or storage in many research area.

Materials


Tube: PP (Polypropylene)/PS (Polystyrene), Compliance with USP Class VI Color: Clear
Cap: HDPE (High Density Polyethylene), Compliance with USP Class VI Color: Green
Sealing ring: TPE (Thermoplastic Elastomer) Color: Gray
Rack: Paper Color: White

Dimensions (Unit:mm)



Features

- Easy-to-read black graduations in ±2% increments
- Contains a large, white frosted writing area
- Autoclavable at 121°C and freezable to -80°C
- Meet the requirements of BSE/TSE
- Latex Statement: This product is latex free
- Non-pyrogenic: The acceptance level for product is 0.05 EU/ml. (Tachypleus Amebocyte Lysate (TAL) Clot-gel test)
- Meet the requirements of RoHS
- Leak-proof
- DNase/RNase free
- Sterilized by gamma irradiation SAL10^{-6}(ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® 15mL Ultra High Performance Centrifuge Tubes

VWR EU EAN. No.	VWR NA Cat. No.	Capacity (mL)	Maximum RCF	Tube Material	Tube Shape	Tube Color	Cap Shape	Cap Color	Sterilization	Packaging Configuration
525-1084	21008-103	15mL	17,000g	PP	Conical-Bottom	Clear	Flat cap with sealing ring	Green	Yes	25/RK, 2 RK/PK, 10 PK/CS, 500/CS
525-1081	21008-089	15mL	17,000g	PP	Conical-Bottom	Clear	Flat cap with sealing ring	Green	No	25/RK, 2 RK/PK, 10 PK/CS, 500/CS
525-1082	21008-214	15mL	17,000g	PP	Conical-Bottom	Clear	Flat cap with sealing ring	Green	No	50/Bag, 10 Bags/CS, 500/CS
525-1085	21008-216	15mL	17,000g	PP	Conical-Bottom	Clear	Flat cap with sealing ring	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1083	21008-105	15mL	17,000g	PP	Conical-Bottom	Clear	Flat cap with sealing ring	Green	No	Tube:500/Bag, Cap:500/Bag, 500/CS
525-1089	93000-026	15mL	17,000g	PP	Conical-Bottom	Clear	Plug cap	Green	Yes	25/RK, 2 RK/PK, 10 PK/CS, 500/CS
525-1086	93000-024	15mL	17,000g	PP	Conical-Bottom	Clear	Plug cap	Green	No	25/RK, 2 RK/PK, 10 PK/CS, 500/CS
525-1090	93000-022	15mL	17,000g	PP	Conical-Bottom	Clear	Plug cap	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1087	93000-020	15mL	17,000g	PP	Conical-Bottom	Clear	Plug cap	Green	No	50/Bag, 10 Bags/CS, 500/CS
525-1088	93000-028	15mL	17,000g	PP	Conical-Bottom	Clear	Plug cap	Green	No	Tube:500/Bag, Cap:500/Bag, 500/CS
525-1091		15mL	20,000g	PP	Conical-Bottom	Clear	Flat cap with sealing ring	Green	Yes	25/RK, 2 RK/PK, 10 PK/CS, 500/CS
525-1092		15mL	20,000g	PP	Conical-Bottom	Clear	Flat cap with sealing ring	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1437		15mL	20,000g	PP	Conical-Bottom	Clear	Flat cap with sealing ring	Green	Yes	Individually wrapped, 100/CS

Technical Data Sheet



VWR® 15mL Ultra High Performance Centrifuge Tubes

Description: 15mL Ultra High Performance Centrifuge Tubes, Flat cap with sealing ring, PS Sterilized/Non-Sterilized

Purpose: Centrifuge tubes are popular in samples centrifuging or storage in many research area.

Materials


Tube: PS (Polystyrene), Compliance with USP Class VI Color: Clear


Cap: HDPE (High Density Polyethylene), Compliance with USP Class VI Color: Green

Sealing ring: TPE (Thermoplastic Elastomer) Color: Gray

Rack: Paper Color: White

Dimensions (Unit:mm)





Features

- Molded graduations stand out distinctly on these optically clear polystyrene tubes
- Colorless and transparent are perfect for biological reaction experiments
- Meet the requirements of BSE/TSE
- Latex Statement: This product is latex free
- Non-pyrogenic: The acceptance level for product is 0.05 EU/ml. (Tachypleus Amebocyte Lysate (TAL) Clot-gel test)
- Meet the requirements of RoHS
- Leak-proof
- DNase/RNase free
- Sterilized by gamma irradiation SAL10^{-6}(ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 15mL Ultra High Performance Centrifuge Tubes

VWR EU EAN. No.	VWR NA Cat. No.	Capacity (ml)	Maximum RCF	Tube Material	Tube Shape	Tube Color	Cap Shape	Cap Color	Sterilization	Packaging Configuration
525-1093	21008-202	15mL	3,000g	PS	Conical-Bottom	Clear	Flat cap with sealing ring	Green	Yes	25/RK, 2 RK/PK, 10 PK/CS, 500/CS
525-1094	21008-212	15mL	3,000g	PS	Conical-Bottom	Clear	Flat cap with sealing ring	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1095	21008-189	15mL	3,000g	PS	Conical-Bottom	Clear	Flat cap with sealing ring	Green	No	Bulk, 1000/CS

Technical Data Sheet



VWR® 50mL Ultra High Performance Centrifuge Tubes

Description: 50mL Ultra High Performance Centrifuge Tubes with Flat or Plug Caps, Flat cap with sealing ring, Sterilized/Non-Sterilized

Purpose: Centrifuge tubes are popular in samples centrifuging or storage in many research area.

Materials


Tube: PP (Polypropylene), Compliance with USP Class VI Color: Clear

Cap: HDPE (High Density Polyethylene), Compliance with USP Class VI Color: Green

Sealing ring: TPE (Thermoplastic Elastomer) Color: Gray

Rack: Paper Color: White

Dimensions (Unit:mm)



Features

- Easy-to-read black graduations in ±2% increments
- Contains a large, white frosted writing area
- Autoclavable at 121℃ and freezable to -80℃
- Meet the requirements of BSE/TSE
- Latex Statement: This product is latex free
- Non-pyrogenic: The acceptance level for product is 0.05 EU/ml. (Tachypleus Amebocyte Lysate (TAL) Clot-gel test)
- Meet the requirements of RoHS
- Leak-proof
- DNase/RNase free
- Sterilized by gamma irradiation SAL10^{-6}(ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system

②

STERILE	R

VWR® 50mL Ultra High Performance Centrifuge Tubes

VWR EU EAN. No.	VWR NA Cat. No.	Capacity (mL)	Maximum RCF	Tube Shape	Tube Color	Cap Shape	Cap Color	Sterilization	Packaging Configuration
525-1098	21008-169	50mL	20,000g	Conical-Bottom	Clear	Flat cap with sealing ring	Green	No	25/RK/PK, 20 PK/CS, 500/CS
525-1107	21008-178	50mL	20,000g	Conical-Bottom	Clear	Flat cap with sealing ring	Green	Yes	25/RK/PK, 20 PK/CS, 500/CS
525-1099	21008-240	50mL	20,000g	Conical-Bottom	Clear	Flat cap with sealing ring	Green	No	50/Bag, 10 Bags/CS, 500/CS
525-1109	21008-242	50mL	20,000g	Conical-Bottom	Clear	Flat cap with sealing ring	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1100	21008-177	50mL	20,000g	Conical-Bottom	Clear	Flat cap with sealing ring	Green	No	Tube:500/Bag, Cap:500/Bag, 500/CS
	37001-468	50mL	20,000g	Conical-Bottom	Clear			No	Tube:500/Bag/CS
525-1111	93000-036	50mL	20,000g	Conical-Bottom	Clear	Plug cap	Green	Yes	25/RK/PK, 20 PK/CS, 500/CS
525-1101	93000-034	50mL	20,000g	Conical-Bottom	Clear	Plug cap	Green	No	25/RK/PK, 20 PK/CS, 500/CS
525-1113	93000-032	50mL	20,000g	Conical-Bottom	Clear	Plug cap	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1103	93000-030	50mL	20,000g	Conical-Bottom	Clear	Plug cap	Green	No	50/Bag, 10 Bags/CS, 500/CS
525-1105	93000-038	50mL	20,000g	Conical-Bottom	Clear	Plug cap	Green	No	Tube:500/Bag, Cap:500/Bag, 500/CS
	82018-048	50mL	10,000g	Freestanding	Clear	Flat cap with sealing ring	Green	No	50/Bag, 10 Bags/CS, 500/CS
525-1097	82018-050	50mL	10,000g	Freestanding	Clear	Flat cap with sealing ring	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1096	82018-052	50mL	10,000g	Freestanding	Clear	Flat cap with sealing ring	Green	No	Tube:500/Bag, Cap:500/Bag, 500/CS
525-1311		50mL	20,000g	Conical-Bottom	Clear	Flat cap with sealing ring	Green	Yes	Individually wrapped, 25/PK, 500/CS

Technical Data Sheet



VWR® Ultra High Performance Light Sensitive Centrifuge Tubes

Description: Ultra High Performance Light Sensitive Centrifuge Tubes with Plug Caps, Sterilized/Non-Sterilized

Purpose: Centrifuge tubes are popular in samples centrifuging or storage in many research area.


Materials


Tube: PP (Polypropylene), Compliance with USP Class VI Color: Amber

Cap: HDPE (High Density Polyethylene), Compliance with USP Class VI Color: Green

Rack: Paper Color: White


Dimensions (Unit:mm)





Features

- Easy-to-read black graduations in ±2% increments
- Contains a large, white frosted writing area
- Autoclavable at 121℃ and freezable to -80℃
- Meet the requirements of BSE/TSE
- Latex Statement: This product is latex free
- Non-pyrogenic: The acceptance level for product is 0.05 EU/ml. (Tachypleus Amebocyte Lysate (TAL) Clot-gel test)
- Meet the requirements of RoHS
- Leak-proof
- Block UV Lighting
- DNase/RNase free
- Sterilized by gamma irradiation SAL10^{-6}(ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® Ultra High Performance Light Sensitive Centrifuge Tubes

VWR EU EAN. No.	VWR NA Cat. No.	Capacity (mL)	Maximum RCF	Tube Shape	Tube Color	Cap Shape	Cap Color	Sterilization	Packaging Configuration
525-1116	89079-543	15mL	12,500g	Conical-Bottom	Amber	Plug cap	Green	Yes	25/RK, 2 RK/PK, 10 PK/CS, 500/CS
525-1117	89079-541	15mL	12,500g	Conical-Bottom	Amber	Plug cap	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1115	89079-538	15mL	12,500g	Conical-Bottom	Amber	Plug cap	Green	No	Tube:500/Bag, Cap:500/Bag, 500/CS
525-1119	89079-537	50mL	12,500g	Conical-Bottom	Amber	Plug cap	Green	Yes	25/RK/PK, 20 PK/CS, 500/CS
525-1120	89079-535	50mL	12,500g	Conical-Bottom	Amber	Plug cap	Green	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1118	89079-532	50mL	12,500g	Conical-Bottom	Amber	Plug cap	Green	No	Tube:500/Bag, Cap:500/Bag, 500/CS

Technical Data Sheet



VWR® Metal-Free Centrifuge Tubes

Description: Metal-Free Centrifuge Tubes with Flat Caps, Sterilized

Purpose: Centrifuge tubes are popular in samples centrifuging or storage in many research area.


Materials

Tube: PP (Polypropylene), Compliance with USP Class VI Color: Clear

Cap: HDPE (High Density Polyethylene), Compliance with USP Class VI Color: Natural

Rack: Paper Color: White

Dimensions (Unit:mm)



Features

- Easy-to-read black graduations in ±2% increments
- Contains a large, white frosted writing area
- Autoclavable at 121℃ and freezable to -80℃
- Meet the requirements of BSE/TSE
- Latex Statement: This product is latex free
- Non-pyrogenic: The acceptance level for product is 0.05 EU/ml. (Tachypleus Amebocyte Lysate (TAL) Clot-gel test)
- Leak-proof
- DNase/RNase free
- Centrifuge Tubes are certified to meet ROHS requirement
- Centrifuge Tubes are certified free of major critical trace metals to refer to the test report on next page
- Sterilized by gamma irradiation SAL10^{-6}(ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system



VWR® Metal-Free Centrifuge Tubes

VWR EU EAN. No.	VWR NA Cat. No.	Capacity (mL)	Maximum RCF	Tube Shape	Tube Color	Cap Shape	Cap Color	Sterilization	Packaging Configuration
525-1121	89049-170	15mL	12,500g	Conical-Bottom	Clear	Flat cap	Natural	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1123	89049-172	15mL	12,500g	Conical-Bottom	Clear	Flat cap	Natural	Yes	25/RK, 2 RK/PK, 10 PK/CS, 500/CS
525-1124	89049-174	50mL	12,500g	Conical-Bottom	Clear	Flat cap	Natural	Yes	50/Bag, 10 Bags/CS, 500/CS
525-1125	89049-176	50mL	12,500g	Conical-Bottom	Clear	Flat cap	Natural	Yes	25/RK/PK, 20 PK/CS, 500/CS

Technical Data Sheet

VWR® Metal-Free Centrifuge Tubes

Test Report

Test items	Measurement unit	Test Results(15mL)	Test Results(50mL)	Test basis
Li	μg/kg	0.0002	0.00007	GB/T30903-2014 ICP-MS
Be	μg/kg	0.0005	0.00002(<0.0001)	
B	μg/kg	0.004	0.019	
Na	μg/kg	0.067	0.060	
Mg	μg/kg	0.029	0.023	
Al	μg/kg	0.24	0.17	
K	μg/kg	0.020	0.018	
Ca	μg/kg	1.9	0.76	
Ti	μg/kg	0.074	0.052	
V	μg/kg	0.0005(<0.0008)	0.0003(<0.0008)	
Cr	μg/kg	0.006	0.009	
Mn	μg/kg	0.008	0.001	
Fe	μg/kg	1.0	0.16	
Co	μg/kg	0.043	0.042	
Ni	μg/kg	0.005	0.005	
Cu	μg/kg	0.046	0.014	
Zn	μg/kg	0.014	0.016	
Ga	μg/kg	0.00003(<0.0001)	0.00002(<0.0001)	
Ge	μg/kg	0.008	0.009	
Sr	μg/kg	0.001	0.001	
Zr	μg/kg	0.0002(<0.0006)	0.00004(<0.0006)	
As	μg/kg	0.001	0.001	
Ag	μg/kg	0.0002	0.0002	
Cd	μg/kg	0.00006(<0.0001)	0.00002(<0.0001)	
Sn	μg/kg	0.0008	0.002	
Sb	μg/kg	0.0007	0.0007	
Ba	μg/kg	0.006	0.0007	
Ta	μg/kg	0.00004(<0.0002)	0.00007(<0.0002)	
W	μg/kg	0.0001(<0.0006)	0.0001(<0.0006)	
Au	μg/kg	0.0001	0.00004(<0.0001)	
Tl	μg/kg	0.0003	0.0003	
Pb	μg/kg	0.003	0.001	
Bi	μg/kg	0.00002(<0.0007)	0.00001(<0.0007)	
Se	μg/kg	0.11	0.073	
Hg	μg/kg	0.004	0.003	
Si	μg/kg	0.20	0.26	
P	μg/kg	0.71	0.57	

Remark:

Preparation and test method of the extracted solution of centrifuge tube: First add purified 5% HF solution into the tube, then tighten the cap and put it in the shaker for extraction in 120rpm for 24 hours. The extracted solution is defined as the test solution which is determined the contents of 37 kinds of reported elements by using ICP-MS.

Technical Data Sheet



VWR® 15mL Polystyrene(PS) Centrifuge Tubes

Description: 15mL Polystyrene(ps) Centrifuge Tubes with Flat Caps, Sterilized/Non-Sterilized

Purpose: Centrifuge tubes are popular in samples centrifuging or storage in many research area.


Materials


Tube: PS (Polystyrene), Compliance with USP Class VI Color: Clear

Cap: HDPE (High Density Polyethylene), Compliance with USP Class VI Color: Green

Gasket: PE (Polyethylene) Color: White

Dimensions (Unit:mm)





Features

- Molded graduations stand out distinctly on these optically clear polystyrene tubes
- Colorless and transparent are perfect for biological reaction experiments
- Meet the requirements of BSE/TSE
- Latex Statement: This product is latex free
- Non-pyrogenic: The acceptance level for product is 0.05 EU/ml. (Tachypleus Amebocyte Lysate (TAL) Clot-gel test)
- Meet the requirements of RoHS
- Leak-proof
- DNase/RNase free
- Sterilized by gamma irradiation SAL 10^{-6} (ISO11137)
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system

VWR® 15mL Polystyrene(PS) Centrifuge Tubes

VWR EU EAN. No.	VWR NA Cat. No.	Capacity (mL)	Maximum RCF	Tube Shape	Tube Color	Cap Shape	Cap Color	Sterilization	Packaging Configuration
	21008-187	15mL	3,000g	Conical-Bottom	Clear			No	Bulk, 1000/CS
	21008-197	15mL	3,000g	Conical-Bottom	Clear	Flat cap	Green	Yes	125/Bag, 4 Bags/CS, 500/CS