

T-Port Extension

With injection port and side clamp

- Small bore infusion extension with side clamp and injection site
- Outer ø: 2.4 mm, inner ø: 1.3 mm, Volume 0.1 ml, Length 10 cm

DESCRIPTION	PRIMING VOLUME (ml)	LENGTH	INNER Ø	UNITS / UOM	CODE#
Small bore ext. with T-port and luer slip	0.1	10 cm	1.3 mm	100	471950
Small bore ext. with T-port and spin lock	0.1	10 cm	1.3 mm	100	471954

LS-Connector

Pressure resistant up to 2 bar

DESCRIPTION		PRODUCT DESCRIPTION	LENGTH	UNITS / UOM	CODE#
LS-2-connector		With flexible tube parts	15 cm	100	4097122
		With back-check valves and clamps	15 cm	50	4093038
LS-4-connector		With flexible tube parts and roller clamp	35 cm	50	4097149

Special Devices

Dosifix®, Twin Spike

SPECIAL DEVICES		MATERIAL	TUBE LENGTH	TUBE LENGTH	INJ. PORT	CODE #
Dosifix® Luer Lock Burette System		PVC,DEHP free	160 cm	63"	I.S.	4037014
Twin Spike Vet Infusion Set for 2 bags		PVC,DEHP free	250 cm	98"	-	4087193