

REV	DESCRIPTION OF CHANGE	DRAWN BY	DATE	ECN
B	REVISED DRAWING	TRJ	08/11/09	20090803102517 58
–	REDRAWN ON DWK TEMPLATE AND STDS	EL	07/19/19	N/A



ACTUAL SIZE



ITEM	QTY.	MATERIAL: DESCRIPTION AND SPECIFICATIONS
1	1	CLEAR COATED OR LACQUERED 3003–H14 ALUMINUM, .008 ±.001 THICK

THIRD ANGLE PROJECTION

TITLE: 13 MM SEAL, CTR TEAR–OUT, ALUM NAT, UNLN

THIS CONFIDENTIAL AND PROPRIETARY DRAWING AND THE INFORMATION IN IT ARE THE EXCLUSIVE PROPERTY OF DWK LIFE SCIENCES. IT IS PROVIDED WITH THE UNDERSTANDING AND INTENTION THAT NOTHING CONTAINED HEREIN WILL BE REPRODUCED OR DISCLOSED TO THIRD PARTIES WITHOUT THE PRIOR EXPRESS WRITTEN CONSENT OF DWK LIFE SCIENCES.

ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF DWK LIFE SCIENCES ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

TOLERANCES
DIMENSIONS IN INCHES UNLESS SPECIFIED
.XX ± .01"
.XXX ± .005"
.XXXX ± .0005"
FRACTIONAL ± 1/64"

DURAN WHEATON KIMBLE

ITEM NUMBER: 224182–01	DRAWING NUMBER: CSD–224182–01	DRAWN BY: TRJ	DRAWN DATE: 10/23/97
SAP NUMBER: 224182–01	OLD MATERIAL NUMBER: N/A	SHEET SIZE: A	SHEET NUMBER: 1 OF 1 / SHEET SCALE: 3:1