

REV	DESCRIPTION OF CHANGE	DRAWN BY	DATE	ECN
-	INITIAL ISSUE	ESM	09/03/04	N/A
-	REDRAW ON DWK TEMPLATE AND STDS	EL	12/01/21	N/A





FINISH DETAIL
SCALE 4:1

NOTES:
1. TUBE DOES NOT REQUIRE THE USP<660> COMPLIANCE.

ITEM	QTY.	MATERIAL: DESCRIPTION AND SPECIFICATIONS
1	1	CLEAR BOROSILICATE GLASS, 33 EXPANSION AS PER ASTM E438 - SEE NOTE 1

THIRD ANGLE PROJECTION

TITLE:

TUBE,CULTURE,PTFE LINED,13-415,13X100M

THIS CONFIDENTIAL AND PROPRIETARY DRAWING AND THE INFORMATION IN IT ARE THE EXCLUSIVE PROPERTY OF DWK LIFE SCIENCES. IT IS PROVIDED WITH THE UNDERSTANDING AND INTENTION THAT NOTHING CONTAINED HEREIN WILL BE REPRODUCED OR DISCLOSED TO THIRD PARTIES WITHOUT THE PRIOR EXPRESS WRITTEN CONSENT OF DWK LIFE SCIENCES.

ELECTRONIC OR PAPER COPIES OF THIS PRINT OUTSIDE OF DWK LIFE SCIENCES ARE UNCONTROLLED AND ARE FOR REFERENCE ONLY UNLESS OTHERWISE STAMPED.

DURAN WHEATON KIMBLE

ITEM NUMBER:	DRAWING NUMBER:	DRAWN BY:	DRAWN DATE:
45066A-13100	(45066A-13100)	ESM	

SAP NUMBER:	OLD MATERIAL NUMBER:	SHEET SIZE:	SHEET NUMBER:	SHEET SCALE:
	N/A	A	1 OF 1	1:1