# Product Data Sheet

## Commercial Grade Glass Hydrometer
## Precision line
Reference 3011FC120/20-qp

For general purpose


Density hydrometer (g/ml)

With conical bottom

Reading below meniscus

Not available for Calibration Certificate






## Specifications

- Range 1200-1300 g/ml
- Division: 0.001 g/ml
- Accuracy: ± 0.001 g/ml
- Calibration temperature : 20°C
- Surface tension : 55 mN/m
- ± 285 mm
- In individual box

Alla France
Z.I. du Bompas - 13, rue du Bompas - Chemillé - 49120 Chemillé-en-Anjou - FRANCE
Tél. + 33 (0)2 41 30 55 08 - Fax +33 (0)2 41 30 34 67
www.allafrance.com