
BRAND GMBH + CO KG

7013 40 96-well Deep-well Plates

Volume:	1.1 ml
Material:	PP
Color:	Colorless
Bottom form:	U-bottom
Surface:	Non-treated, non-sterile
Storage:	In closed original cartons, dry, at room temperature

Dimensions:



L	127.79 mm
W	85.49 mm
H	27.20 mm
Hs	2.50 mm
X1	12.40 mm
Xn	113.38 mm

Y1	11.24 mm
Yn	74.24 mm
P	9.00 mm
C	12
R	8