Laboratory beakers/Griffin beakers, PP



Laboratory beakers according to ISO 7056, based on DIN EN ISO 3819 (glass, formerly DIN 12331) for various uses.
When bought in bulk, individual print is possible.

- PP, highly transparent
- Blue graduation printed on
- Shatterproof
- Autoclavable

Laboratory beakers/Griffin beakers, PP

Capacity ml	Division ml	Scaling from ml	Dia. mm	Height mm	Price €	Quantity scale from	Scale price €	Item no.	
25	1	5	34	49		12		7122-0025	N3
50	5	10	43	60		12		7122-0050	N3
100	5	10	51	72		12		7122-0100	N3
150	5	10	58	82		12		7122-0150	N3
250	5	20	70	95		6		7122-0250	N3
400	10	20	80	110		6		7122-0400	N3
500	50	50	94	118		6		7122-0500	N3
600	25	50	93	125		12		7122-0600	N3
1000	50	100	108	148		6		7122-1000	N3
2000	100	200	140	188		6		7122-2000	N3
3000	250	500	170	212		4		7122-3000	N3
5000	500	500	203	247		4		7122-5000	N3

Graduated beakers, stackable



Graduated beaker with graduation printed in blue, open handle, conical. Beakers can be stacked and therefore save space during storage.

- PP, transparent
- Sterilisable
- Shatterproof
- Good chemical resistance

Graduated beakers, stackable

Capacity ml	Division ml	Height mm	Dia. top mm	Price €	Quantity scale from	Scale price €	Item no.	
1000	10	167	125		12		7110-1000	N3
2000	20	212	148		12		7110-2000	N3