

Vitrex Steel

General Purpose Capillary Blood Lancet

Vitrex Steel lancets are traditional handheld, flat and stainless steel blood lancets for general purpose use. The steel lancets are available in standard and in allergy version.

LOW COST
Simple, flat, steel lancet.

SAMPLE VOLUME
Excellent blood flow.

HIGH QUALITY
Hardened surgical steel.

MINIMUM DISCOMFORT
Modern assymetric design.

STERILE
Sterilised by gamma radiation to EN 556 and EN ISO 11137.

Technical Specifications

Product	Description	Article No.
Vitrex Steel - standard	Flat blade sterile lancet	350213
Vitrex Steel - allergy	Flat blade sterile lancet	365213

Shelf life 5 years



Transit packaging Pallet

Packaging:

Primary packaging:	200 pcs.
Secondary packaging:	10,000 pcs.
Transit packaging:	40,000 pcs.
Pallet:	720,000 pcs.