Technical Data Sheet

VWR® PET BOTTLES



Description: VWR® Square media 1000ML(33.8OZ) bottle, white cap, Sterile, Tray package Bottles are packaged individually wrapped within two outer packs per case

Purpose: To use on packing, storage or transport for powder and liquid

Bottle Material: Polyethylene terephthalate (PET), Color: Transparent

Cap Material: Polyethylene (PE), Color: White

Features

- EB Sterilization SAL:10^{-6} (ISO11137)
- Sterile Package, light, solid and good sealing feature
- Products are manufactured in class 100000 cleaning room
- Shelf Life: 3 years after month of production.
- Products meet the requirements of BSE/TSE
- Products meet the BPA Statement
- Products meet the Heavy metals - ROHS
- Engraved graduation
- Non-pyrogenic
- Manufactured under EN ISO 13485 quality management system

VWR® PET BOTTLES

NA Part No.	EU Part No	Size(mm)			Neck I.D.（mm）	Cap Size（mm）	Outer Diameter (mm)	EB Sterile	Individually Pack	Packing
		Thickness	Width	Height				SAL 10^{-6}		
76302-840-Jar	VWRI391-0631-Jar	0.7±0.5	93.0±1.5	212.3±1.5	29.0±0.3	37.8±0.3	93.0±1.5	Yes	Yes	24pcs/case (12pcs/tray,2trays/case)
		Width		Height						
76302-840-Lid	VWRI391-0631-Lid	41.9±0.3		29.6±0.4						

Drawing of bottle

