Product name : **CytoSep**	Catalogue no.: **M964-10FW**
Edition date : 2010-11-05 Approved by : Annette Roy, QA coordinator	Revision : 2014-06-19
This document replaces any previous version	

↺ Product identification

- Designation: **Disposable cytology funnel with white filter paper**

↺ Production Information

- Producer's name: **Simport *Scientific***
- Address: 2588 Bernard Pilon, Beloeil (Québec), J3G 4S5
- Tel. no / web site: (450) 464-1723 / www.Simport.com

↺ Packaging

- Case: Cardboard box containing 10 packs of 50 units, for a total of 500 units.

↺ Traceability.

- Lot no
  - ↪ Composition: 8 digits
  - ↪ Location: (The lot number can be found in different locations)
    ① On exterior case label.
    ② On inserted label inside case.

↺ Standard conformity.

- **MIL-STD-105E**: Sampling and inspection procedures
- **Coneg**: Plastic resin is in conformity with CONEG standards for heavy metals
- **FDA**: Plastic resin conforms to 21CFR 177.1520
  Filter paper conforms to 21CFR 1476.170 and 176.180

Product name : **CytoSep**	Catalogue no.: **M964-10FW**

↺Quality assurance

- Clear, no presence of contamination in plastic
- Visual attributes

↺ Specifications

- Single cytology funnel for the Shandon Cytospin Cytocentrifuges
- Material: Polypropylene funnel with cap and white filter paper

↺ Storage conditions

- Store in cool dry place.
- Avoid temperature variations and humidity.
- Protect from any possible contamination.

↺ Recommended use

- Careful for spills when using biohazard material and/or chemical reagents.
- Follow chemical resistance chart recommendations

OWNERSHIP OF MATERIALS:
Materials and information contained in this document are copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.