Chemical spoons



Chemical spoon with spatula handle, stainless steel

with spatula handle

Material	A (Length)	a (Length)	B (Width)	D/Ø	F x E	Article no.
18/10 Stainless	120	30	20	1	31x23	3420
18/10 Stainless	150	32	22	1	38x27	3421
18/10 Stainless	180	32	22	1	40x28	3422
18/10 Stainless	210	32	22	1	40x28	3423
18/10 Stainless	250	45	32	2	50x35	3424
18/10 Stainless	300	50	37	2	55x41	3425
Ni 99,5%	120	30	20	1	31x23	3430
Ni 99,5%	150	32	22	1	38x27	3431
Ni 99,5%	180	32	22	1	40x28	3432
Ni 99,5%	210	32	22	1	40x28	3433
PTFE coating	180	32	22	1	40x28	3750
Titanium	180	32	22	1	40x28	3920
Carbon coating	180	32	22	1	40x28	3426

Bochem Instrumente GmbH · Industriestraße 3 · D-35781 Weilburg
Tel. +49 (0) 6471 9282-0 · Fax +49 (0) 6471 9282-30 · info@bochem.de · www.bochem.com