BAKERBOND spe, Silica Gel

BAKERBOND spe

Disposable Extraction Columns; Glass

Material No: 7337-07
Effective from: 2019/11/28

# Product Specification

| Test | Specification |
|---|---|
| Average Particle Size, um | about 50 |
| Average Pore Diameter, A | about 60 |
| Loss on drying (at 180°C) | <= 6.0 % |
| Particle size distribution | 35 - 70 um |
| pH of 10% aquous susp. at 25°C | 6 - 8 |
| The following results are based on the silica gel raw material: | |
| Average Particle Diameter, um (APD) | 47 - 60 |
| Specific Surface Area, m2/g | For Information Only |
| Particle size d90 [um] | <= 105.0 |
| Particle size d10 [um] | >= 33.0 |

For Laboratory, Research or Manufacturing use.
8mL Solid Phase Extraction Columns
1000 mg per Column

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000