boettger

Technical Data (TD)

Faeces Container 15 ml, Screw Closure					
REF	07-023-2007				

1.	**Description / Specification**	**Tube**	**Cap**	**Spoon**
1.01	Description	Faeces Container 20/76, Screw Cap w/Spoon		
1.02	Dimension	See *Customer Drawing*		
1.03	Volume	15 ml		
1.04	Material / Resin	Polypropylene (PP)	Polyethylene (HDPE)	Polystyrene (HIPS)
	Grade of Resin	MAT120	MAT131	MAT103
	Properties of Resin	See *www.boettger.com*		
1.05	Colour	Translucent	Brown	Opaque
	Grade of Pigment		PIG117	
1.06	Hygiene Level	Standard hygienic		
1.07	Quality Control	Testing of attributive and variable characteristics in accordance with the valid specification		
1.08	Usage / Purpose	General laboratory product, Collection and evaluation of stool samples		
1.09	Conformities / Regulations	See *Classification*		
1.10	Product Stewardship	See *Product Stewardship* of Grades of Resins/Pigments		
1.11	Storage	in closed original packaging, dry, at room temperature		
1.12	Transport	Dry, at room temperature, handle with care, do not throw or push cartons, do not stack pallets		
1.13	Other information	For single use only		

2.	**Features**	
2.01	Basic features	-
2.02	Temperature range	-20°C to 60°C
2.03	Autoclavability	No
2.04	Centrifugation	N/A
2.05	Chemical Resistance	See *www.boettger.com*
2.06	Expiry	5 years from LOT (YYYY-MM)
2.07	Other information	Leakproof corresponding to EN ISO 14254 D

3.	**Packaging**			
3.01	Pcs/Bag	250		
3.02	Pcs/Ctn	1000		
3.03	Pcs/Pallet	24000		
3.04	LOT / VP	WWYY (week, year of production) / XX (packaging)		
3.05	Other information	See *Measures*		

4.	**Other Information**	
4.01	Codes and Numbers	See *Codes and Numbers*
4.02	Symbols	See *Symbols*
4.03	CH-REP	Swiss AR Services AG, Industriestr. 47, CH-6300 Zug CHRN-AR-20000807

Confidential: Information contained in this document or drawing is confidential and proprietory to Paul Boettger GmbH & Co. KG. This document may not be reproduced for any reason without written permission from Paul Boettger GmbH & Co. KG. All rights of design, invention and copyright are reserved.

5.	Label (Example) according to IVDR 2017/746/EC



6.	Printing on Subpackaging (Example)



All informations without guarantee. Above mentioned products are manufactured according to company specifications and our certified Quality Management System corresponding to DIN EN ISO 13485. The quality check is carried out on a sample basis; an individual examination is not possible. Deficiencies parts may occur. We accept no liability for consequential damages.