APPLICATIONS:
For tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in an aseptic environment.

FEATURES:
- Sterility assurance on individually wrapped pipet
  - Contamination control
- Consistent performance
  - Reliability, contamination control
- Aseptic presentation
  - Contamination control
- Non-pyrogenic
  - Low endotoxins
- Crystal grade polystyrene
  - Can clearly see fluid sample
- Non-wettable
  - Minimum meniscus
- Individually wrapped packaging: Cellulose-based, direct-seal paper on one side and the bottom web is polyethylene plastic material.

REGULATORY COMPLIANCE:
Falcon® products are manufactured under the current ISO 13485 Standards, and the current FDA Quality System Regulation 21 CFR Section 820.

STERILITY:
Product has been sterilized and dosimetrically released per the requirements of ANSI/AAMI/ISO 11137, "Sterilization of Health Care Products – Radiation." Product meets a minimum Sterility Assurance Level (SAL) of 10^{-6}.

Warranty/Disclaimer: Unless otherwise specified, all products are For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.


For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.
© 2022 Corning Incorporated. All rights reserved.

Falcon® 2 mL Aspirating Pipet

Material: Polystyrene

Related Catalog Numbers:	Cat. No.	Packout Description
Sterile	357558	Individually wrapped, 50/ Pk 200/case

DIMENSIONS = XXX mm (XXX inches)	
TIP: Length: = N/A Interior Diameter: = 1.625 (0.64) Outside Diameter: = 3.73 (.147)	**PIPET SIZE:** 2 mL **INCREMENTS:** N/A
MOUTHPIECE: Length: = N/A Interior Diameter: = 4.165 (.164) Outside Diameter: = 6.27 (.247) Plug Length: = N/A	**GRADUATIONS:** N/A **OVERFILL CAPACITY:** N/A **COLOR CODE:** N/A
EXTRUSION: Interior Diameter: = N/A Outside Diameter: = N/A Length: = N/A	**GENERAL DIMENSIONS:** A = 282.58 (11.125) B = 3.73 (.147) Diameter C = 6.27 (.247) Diameter D = 6.27 (.247) Diameter



Warranty/Disclaimer: Unless otherwise specified, all products are For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.
© 2022 Corning Incorporated. All rights reserved.