Technical sheet

Kit Myeloperoxydase

Product code: 361610-0000

Composition:

• **Alphanaphtol Solution**	code: 320750-0030	1 x 30 mL
• **Formaldehyde-Ethanol solution**	code: 310850-0125	1 x 125 mL
• **Mayer Haematoxylin**	code: 361620-0125	1 x 125 mL
• **Hydrogen peroxide**	code: 300980-0010	1 x 10 mL
• **Pyronin Solution**	code: 361800-0125	1 x 125 mL

Applications:

Hematology: Kit Myeloperoxidase allows a staining with Alphanaphtol-Pyronin to points out peroxidase granulations.

Peroxidase granulations form an insoluble and bright red complex with Pyronin in presence of Alphanaphtol and Hydrogen Peroxide.

Mayer Haematoxylin stains the leukocytes nuclei bleu.

Mise à jour du 18/03/2022