Pressure sprayer







The pressure sprayer features rugged industrial quality for use in the chemical and pharmaceutical industry, for example disinfecting, cleaning, plant protection, thin-layer chromatography, technical and sanitary uses.

The practical 360° spray function also allows for upside-down spraying with the Turn'n'Spray pressure sprayer, making it easy to reach even places that are difficult to access. The sprayer can be used for many purposes, such as disinfecting and cleaning, or also for spraying plant protection products.

- Adjustable spray jet, nozzle diameter 0.8 mm
- Temperature of the liquid: 0 - 30 ° C
- Capacity of 1.5 litres
- Parts in contact with media PE, PP, FKM, AISI 316
- Upside-down spray function

LaboPlast

Pressure sprayer

Capacity ml	Container dia. mm	Height with cap mm	Item no.	
1500	140	305	0309-0100	N4

Pressure sprayer Turn'n'Spray

Capacity ml	Container dia. mm	Height with cap mm	Item no.	
1500	140	305	0309-0120	N4

Accessories

Description	Item no.	
Spraying lance extension for pressure sprayer, 30 cm	0309-0101	N4

Pressure sprayer Food



The food pressure sprayer has a transparent viewing window to make it easy to read the fill level. A rubber-coated ring stand ensures that the sprayer can be set down securely. The built-in safety valve protects against overpressure (max. 3 bar). The adjustable nozzle nut makes it possible to regulate the spray jet individually – from a fine spray mist to a powerful jet. The ergonomic handle provides a comfortable hold on the pressure sprayer. The food-safe spray bottle is specially designed for use in the food industry and is well suited for spraying different media including oils, greases or flavouring substances. It can also be used to easily atomise cleaning and disinfecting agents.

- Parts in contact with media PE, PP, VMQ, AISI 304
- Capacity of 1.2 litres

Pressure sprayer Food

Capacity ml	Container dia. mm	Height with cap mm	Item no.	
1200	117	315	0309-0130	N4