Artikel-Nr.: 1-0600-0250-22

Hexagonal "Meplat" glass bottle; clear glass; thread GL 22; vol. 250 ml



Dimensions	
Volume	250 ml
Height	183 mm
weight	255 g
Thread type	GL 22

Material / Propertis	
Quantity	36, 1152
Materials	Glass - Hydrolytische Klasse 3 (Typ III)
Colour	clear
Product Carbon Footprint	591.3 g CO2e / St.
Warning sign	No
Visual strip	No

Accessories

Item number	Item name
8-1001-22	Srew cap closure; PP black; with LDPE insert for thread GL 22;
8-1002-22-W	Screw cap closure; PPN white; with foam inlay for thread GL 22;
8-1004-22	screw cap closure; PP black; with PTFE-inlay 2mm; for thread GL22;

Otto-Hahn-Straße 19, D-68169 Mannheim | Telefon: +49 (0)621-322 72-0
info@rixius.com | www.rixius.com | www.berlinpackaging.eu