# Technical Data Sheet

Product Name: **Cytoslides two circles uncoated**
Article Number: **309-100-0**

This document replaces any previous version.

**Product Description**

- Single use Cytoslides with two circles and uncoated

**Packaging**

- Case: 1 box of 100 units

**Product Specifications**

Compatible with all **Tharmac Cellspins, Shandon Cytospin® Cytocentrifuge** and Easy Rotor

**Materials**

- White glass

**Standards and Conformity**

- CE: The product bears the CE mark

**Quality Assurance**

- Clear, no presence of contamination on the slide
- Visual attributes
- Dimension measurements
- Assembly fit

**Traceability**

Lot No. Composition: 5 digits

- The Lot number can be found on label of each package with 100 units.

**Storage Conditions**

- Store in a warehouse at room temperature
- Avoid humidity and temperature variations
- Protect from any possible contamination

**Recommended Use**

Use in **Tharmac Cellspins, Shandon Cytospin® Cytocentrifuge**, follow the equipment manufacturer's instructions for use.



**THARMAC GmbH**
Im Dachsstück 15 » DE-65549 Limburg an der Lahn » PHONE: +49 (0) 6431 5 84 88-0 » FAX: +49 (0) 6431 5 84 88-28
MAIL: info@tharmac.de » WEB: www.tharmac.de