

thermo scientific

Product Specifications

Product **44565**
Chlorobenzene,99.9%,for HPLC

General Product Data

Version	1
CAS No	108-90-7
Molecular weight	112.56
Molecular formula	C6 H5 Cl
Linear formula	C6H5Cl
Flash point (°C)	23

Product Specifications

Appearance (Form)	Clear liquid
Color scale	=<10 APHA
GC	>=99.85 %
Water	=<0.1 %
UV	at 288 nm A: =<1.0
	at 290 nm A: =<0.40
	at 300 nm A: =<0.05
	at 360 nm A: =<0.01
Residue after evaporation	=<0.0005 %

Product Specifications are subject to amendment and may change over time.

Acros Organics
ENA23, zone 1, nr 1350, Janssen Pharmaceuticalaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: http://www.thermofisher.com
1 Reagent Lane, Fair Lawn, NJ 07410,USA
Fax 201-796-1329