boettger

Technical Data (TD)

Micro Tube PP 1,5 ml w/o Cap					
REF	02-013-0100				

1.	**Description / Specification**			
1.01	Description	Micro Tube 1,5 without cap		
1.02	Dimension	See *Customer Drawing*		
1.03	Volume	1,5 ml		
1.04	Material / Resin Properties of Resin	Polypropylene (RCP)		
		MAT120		
		See *www.boettger.com*		
1.05	Colour Grade of Pigment	**02-013-0100**		
		Translucent		
1.06	Hygiene Level	Standard hygienic		
1.07	Quality Control	Testing of attributive and variable characteristics in accordance with the valid specification		
1.08	Storage	To be stored in closed original packaging: Dry, at room temperature		
1.09	Usage / Purpose	General laboratory product, Collection and evaluation of samples		
1.10	Conformities / Regulations	See *Classification*		
1.11	Product Stewardship	See *Product Stewardship* of Grades of Resins/Pigments		
1.12	Other information	For single use only Also suitable for Eppendorf centrifuges/systems		
2.	**Features**			
2.01	Basic features	-		
2.02	Temperature range	-70°C to 121°C		
2.03	Autoclavability	Yes		
2.04	Centrifugation	18000 g		
2.05	Chemical Resistance	See *www.boettger.com*		
2.06	Expiry	5 years from LOT (YYYY-MM)		
2.07	Other information	-		
3.	**Packaging**			
3.01	Pcs/Bag	1000		
3.02	Pcs/Ctn	10000		
3.03	Pcs/Pallet	240000		
3.04	LOT / VP	WWYY (week, year of production) / XX (packaging)		
3.05	Other information	See *Measures*		
4.	**Other Information**			
4.01	Codes and Numbers	See *Codes and Numbers*		
4.02	Symbols	See *Symbols*		

All informations without guarantee. Above mentioned products are manufactured according to company specifications and our certified Quality Management System corresponding to DIN EN ISO 13485. The quality check is carried out on a sample basis; an individual examination is not possible. Deficiencies parts may occur. We accept no liability for consequential damages.

Confidential: Information contained in this document or drawing is confidential and proprietory to Paul Boettger GmbH & Co. KG. This document may not be reproduced for any reason without written permission from Paul Boettger GmbH & Co. KG. All rights of design, invention and copyright are reserved.