# Forceps curved, blunt points

- curved
- blunt points

---



| Article no. | Material | Length mm | Weight |
|---|---|---|---|
| 3041 | nickel plated steel | 115 | 19 g |
| 3042 | nickel plated steel | 130 | 23 g |
| 3043 | nickel plated steel | 145 | 26 g |
| 3051 | stainless steel 1.4301 | 115 | 19 g |
| 3052 | stainless steel 1.4301 | 130 | 23 g |
| 3053 | stainless steel 1.4301 | 145 | 26 g |

**Carl Friedrich Usbeck KG**
Industriestr. 12-14
42477 Radevormwald
Germany

Phone: +49 (02195) 91 18-0
Fax: +49 (02195) 91 18-40
E-Mail: info@usbeck.eu