Product Specification Sheet

# VWR PCR PLATES, 96-WELL, SKIRTED non-sterile, alpha-numeric labeling, natural-color

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

*VWR International certifie que les produits mentionés respectent toutes les specifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.*

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

VWR Catalog Number : **211-0297**

Product Dimensions: **Length: 128.016mm (5.04")**
**Width: 85.598mm (3.37")**
**Height: 15.494mm (.61")**

Packaging: **10 plates per pack**
**10 packs per case (100 plates total)**

Pyrogens: **Certified testing is performed lot by lot. These tests include Nuclease(Rnase, Dnase) and Endotoxin (pyrogen). Certificates are on file at the production facility.**

Performance Testing: **Each manufacturing lot is sampled and tested in accordance with Standard Operating Procedures.**

***Disclaimer:*** *VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.*