MACHEREY-NAGEL

REF 484000

Filter paper circles, MN 827 ATD, Antibiotic resistance testing







https://qr.mn-net.com/qr/(241)484000

MN 872 ATD products are used in testing the resistance of pathogens to antibiotics. For this test the filter paper can be impregnated with the antibiotic to be tested and placed on the inoculated nutrient medium. Depending on the effectiveness, a smaller or larger zone of inhibition is formed.

Technical data

Parameter*	Value*
Grade	MN 827 ATD
Format	Filter paper circles
Surface	Smooth
Quality	Antibiotic resistance testing
Basis weight	270 g/m^2
Filtration speed description	Fast
Filtration speed	12 s
Thickness	0.68 mm / 26.8 mil
Material	Cellulose
Color	White
Diameter	6 mm
Hazardous material	No Information available

*These data are typical values and must not be regarded as specification

MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com
FR +33 388 68 22 68 sales-fr@mn-net.com
CH +41 62 388 55 00 sales-ch@mn-net.com
US +1 888 321 62 24 sales-us@mn-net.com