Swann-Morton® LIMITED

Penn Works
Owlerton Green
Sheffield, S6 2BJ
Telephone: (0114) 234 4231
Fax: (0114) 231 4966
Email: qa@swann-morton.com
Web: www.swann-morton.com

BLADE REMOVER SPECIFICATIONS

Sterile supplied individually wrapped in boxes of 50
Non-Sterile supplied in boxes of 100 or 500

COMPONENT	MATERIAL
BLADE REMOVER	Phillips K Resin / Grade K R01 butadiene Styrene
TOP WEB - STERILE	Pet / PE MLP12A / 38V (nylon poly)
BOTTOM PAPER - STERILE	60gsm Bleached kraft paper with 10gsm grid laquer
UNIT CONTAINER BOXES STERILE	Coated white lined chipboard with litho machine varnish
UNIT CONTAINER BOXES NON STERILE	Coated white lined chipboard with UV varnish. Non Sterile Blade Removers are double poly bagged.
NON STERILE POLY BAG	Polythene

All components are manufactured in the E.U.

All packs and boxes are marked with the current lot number and expiry dates for traceability.

Sterile products are sterilised by Cobalt 60 Gamma Irradiation by:-Swann-Morton (Services) Limited, Owlerton Green, Sheffield S6 2BJ