

Product Specifications

Product **154630025**
1-Chlorobutane,99+%,pure

General Product Data

Version	3
CAS No	109-69-3
Molecular weight	92.57
Molecular formula	C4 H9 Cl
Linear formula	CH3(CH2)3Cl
Flash point (°C)	-12

Product Specifications

Appearance (Color)	Clear colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms
Refractive index	1.4005 to 1.4035 (20°C, 589 nm)
GC	>=99.0 %

Product Specifications are subject to amendment and may change over time.

Acros Organics
ENA23, zone 1, nr 1350, Janssen Pharmaceuticalaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: http://www.acros.com
1 Reagent Lane, Fair Lawn, NJ 07410,USA
Fax 201-796-1329