BAKERBOND spe, Octadecyl (C18)

BAKERBOND spe

Disposable Extraction Columns; Glass

Material No: 7334-03
Effective from: 2015/06/10

# Product Specification

| Test | Specification |
|---|---|
| The following results are based on the bulk bonded phase: | |
| Appearance | passes test |
| Capacity Factor: Butylbenzene | For Information Only |
| Capacity Factor: Ethylbenzene | For Information Only |
| Capacity Factor: Propylbenzene | For Information Only |
| Carbon (C)(%) | For Information Only |
| Desisopropyl Atrazine (10%), ml | 110 - 150 |
| Surface Coverage (ug/m2), Based on %C | 320 - 350 |
| The following results are based on the silica gel raw material: | |
| Average Particle Diameter, um (APD) - via Malvern (100 mm lens) | 47 - 60 |
| The following results are based on the solid phase extraction columns: | |
| Extraction Residue | <= 0.2 % |

For Laboratory, Research or Manufacturing use.
3mL Solid Phase Extraction Columns
500 mg per Column

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000