Nasco SAMPLING
901 Janesville Avenue, Fort Atkinson, WI 53538-0901

Whirl-Pak® Sample Bag
Material Specification Sheet
B01364WA Stand-Up Bag

Size: 3" W x 7 ¼" L, 7.5 x 18.5 cm.

Capacity: 4 oz, 118 ml.

Film: Blend of linear low and low density virgin polyethylene, 2.5 mil, .064 mm thick, FDA approved for food contact.

Construction: A .75", 1.91 cm bottom gusset makes the bag self standing.

Closing Tab: Two round aluminum wires covered with plastic tape by a patented process that also covers the wire ends, reduces bag puncture, forms leak-proof seal when closed.

Pull Tabs: Made of Tyvek®.

Write-On Strip: Horizontal band of thermoset white ink printed on film, compliant with CONEG resolution for heavy metals and European Union Directive for hazardous substances.

Sterilization: Sterilized with ethylene oxide gas, documentation of sterility is available upon request.

ISO Certification: Manufactured under a quality management system certified to **ISO9001 & ISO13485.**

Packaging: 500 bags/box in a liner bag twist-tied shut.

Condition of Use: Not to be used at temperatures above 180F, 82C.