Technical Data Sheet

Article number	**41325**
Description	**Wide neck bottles, amber glass with thread and screw cap**

Material bottle	Soda-silicate glass, amber
Material cap:	polyurea moulding, black, with lupolen insert, screw caps are heat-resistant up to +110 °C/ +230 °F

Article no.	Capacity	Weight	Height	Diameter	Thread
41325006	25 ml	40 g	57,0 mm	36,0 mm	GL 28
41325012	30 ml	47 g	68,0 mm	36,0 mm	GL 32
41325018	50 ml	65 g	75,0 mm	44,0 mm	GL 32
41325024	75 ml	85 g	79,0 mm	48,5 mm	GL 40
41325030	100 ml	100 g	92,5 mm	50,0 mm	GL 40
41325036	125 ml	110 g	98,0 mm	54,0 mm	GL 40
41325042	150 ml	138 g	104,0 mm	57,0 mm	GL 45
41325048	200 ml	185 g	109,0 mm	63,0 mm	GL 55
41325054	250 ml	218 g	113,0 mm	70,0 mm	GL 55
41325060	300 ml	230 g	127,0 mm	70,0 mm	GL 55
41325066	500 ml	335 g	154,0 mm	84,0 mm	GL 55
41325072	1000 ml	517 g	180,5 mm	103,5 mm	GL 68

Tel.: +49 9779/808-0
Fax: +49 9779/808-88

Stettener Str. 22-24
97647 Sondheim/Rhön · Germany

info@assistent.eu
www.assistent.eu

Glaswarenfabrik Karl Hecht GmbH & Co KG | Amtsgericht Schweinfurt HRA 8063 | USt.IdNr./VAT-No.: DE132198535
Persönlich haftende Gesellschafterin: Hecht GmbH | AG Schweinfurt HRB1160
Geschäftsführer: Karl-Heinz Hecht, Thomas Diepholtz | Einzelprokura: Tonja Müller