Buffer Solution

Baker Analyzed

pH 10.00 [Borate] Blue

Material No: 5658.0500
Effective from: 2016/03/09

Product Specification

Test	Specification
pH at 25°C	9.95 - 10.05

For Laboratory, Research or Manufacturing use.
This buffer is tracable to NIST Standard Reference Buffer.Contains Boric Acid / Potassium Chloride / Sodium Hydroxide.Contains 0.01% Sodium Benzoate.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000