BAKERBOND spe, Amino/Carbon

BAKERBOND spe

Disposable Extraction Columns

Material No: 7450-07
Effective from: 2016/07/22

Product Specification

Test	Specification
Appearance	Small black spheres
Bulk Density (g/ml)	For Information Only
Loss on Drying (%)	For Information Only
Particle size d10 [um]	For Information Only
Particle size d50 [um]	For Information Only
Particle size d90 [um]	For Information Only
Specific Surface Area, by Nitrogen adsorption (BET-method), m2/g	For Information Only
Amino (NH2):	
Carbon (C)(%)	For Information Only
Nitrogen (N)(%)	For Information Only
Extraction Residue	<= 0.2 %
Surface Coverage (ug/m2), Based on %C	For Information Only
Surface Coverage (ug/m2), Based on %N	For Information Only
Particle Size Distribution: APD, microns	47 - 61 microns
Specific Surface Area, m2/g	For Information Only

For Laboratory, Research or Manufacturing use.
6mL Solid Phase Extraction Columns
Per column: 500 mg Amino / 500 mg Carbon

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000