Graduated beakers with handle







Graduated beakers with handle, transparent graduation





Graduated beakers with handle, blue graduation

Graduated beakers made of PP are virtually indestructible and highly transparent.
Graduated beakers are perfectly suited as advertising means for chemical products. The smooth surface enables print of a logo or an advertising slogan.

- PP, transparent
- With indestructible graduation

Graduated beakers with handle, blue graduation

Capacity ml	Dia. top mm	Dia. bottom mm	Height mm	Price €	Quantity scale from	Scale price €	Item no.	
100	50	40	80		24		7104-0100	N3
250	70	50	120		12		7104-0250	N3
500	91	75	140		12		7104-0500	N3
1000	116	95	181		6		7104-1000	N3
2000	150	125	213		6		7104-2000	N3
3000	170	120	242		6		7104-3000	N3
5000	210	160	270		6		7104-5000	N3

Graduated beakers with handle, transparent graduation

Capacity ml	Dia. top mm	Dia. bottom mm	Height mm	Price €	Quantity scale from	Scale price €	Item no.	
100	50	40	80		24		7102-0100	N3
250	70	50	120		12		7102-0250	N3
500	91	75	140		12		7102-0500	N3
1000	116	95	181		6		7102-1000	N3
2000	150	125	213		6		7102-2000	N3
3000	170	120	242		6		7102-3000	N3
5000	210	160	270		6		7102-5000	N3