
3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium pyrophosphate tetrabasic - ≥95%

Product Number:	**P8010**
CAS Number:	7722-88-5
MDL:	MFCD00003513
Formula:	Na4O7P2
Formula Weight:	265.90 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity) 7.6 g plus 100 mL H2O	Clear
% Purity (Titration)	≥ 95
Recommended Retest Period 4 Years	----------

Specification: PRD.0.ZQ5.10000043281

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.