



Canister for storage or transport. Canister with UN transport approval and tamper-evident seal.

- HDPE
- UN approval 3H1/X1.9/250, see page 261
- Screw lid with drain cock available as accessory

Bulk prices

Jerrycans

Capacity l	L x W x H mm	Weight g	Neck dia. inside mm	Thread	Colour	Price €	Quantity scale from	Scale price €	Item no.
5	190 x 145 x 251	250	40,5	DIN50	transparent		16		1427-0005 N3
10	230 x 190 x 310	450	40,5	DIN50	transparent		10		1427-0010 N3
20	290 x 246 x 385	950	48,5	DIN60	blue		6		1427-0020 N3
30	320 x 284 x 435	1250	48,5	DIN60	blue		6		1427-0030 N3

Accessories

Description	Price €	Item no.
Screw closure DIN50 with spigot for jerrycans		1427-9102 N3
Screw closure DIN60 with spigot for jerrycans		1427-9100 N3
Original screw cap DIN50		1427-9105 N3
Original screw cap DIN60		1427-9106 N3