## ViscoDispo







The ViscoDispo disposable sampler can be used to take samples of high-viscosity media up to 100.000 mPas. The ergonomically shaped handle and the very stable drawbar make it easy to draw even viscous substances.

- Ergonomic handle
- No seals
- HDPE
- Clean room manufactured, individually packaged
- Optionally sterilised by gamma rays
- Opening diameter 10 mm
- Ø 21 mm

## PowderDispo



PowderDispo is the ideal disposable sampler for cross-section samples of powdery or slightly cohesive substances. A complete all-layer sample can be taken from bags or smaller containers with the lance.

- HDPE
- Clean room manufactured, individually packaged
- Optionally sterilised by gamma rays
- Ø 21 mm

### ViscoDispo

| Length mm | Volume ml | Model | Pieces per pack | Price € | Quantity scale from | Scale price € | Item no. | |
|---|---|---|---|---|---|---|---|---|
| 500 | 100 | standard | 20 | | 10 | | 5393-3320 | N4 |
| 500 | 100 | sterilised | 20 | | 10 | | 5393-3321 | N4 |
| 1000 | 200 | standard | 20 | | 10 | | 5393-3330 | N4 |
| 1000 | 200 | sterilised | 20 | | 10 | | 5393-3331 | N4 |

### PowderDispo

| Length mm | Volume ml | Model | Pieces per pack | Price € | Quantity scale from | Scale price € | Item no. | |
|---|---|---|---|---|---|---|---|---|
| 300 | 60 | standard | 20 | | 10 | | 5393-5510 | N4 |
| 300 | 60 | sterilised | 20 | | 10 | | 5393-5511 | N4 |
| 500 | 100 | standard | 20 | | 10 | | 5393-5522 | N4 |
| 500 | 100 | sterilised | 20 | | 10 | | 5393-5521 | N4 |