Data Sheet Eye wash bottle BARIKOS KS

BR-Article-No.:	954004
Trade name:	BARIKOS KS
Product-description:	Eye wash bottle according EN 15154-4 Non-perishable water (gamma irradiated)
Kind of product:	First-aid equipment for eye wash
Kind of use:	Eye rinsing while standing / sitting or lying down as a first self-help measure for rinsing out foreign objects and foreign substances
Description:	Transparent plastic bottle with printed on pictogram, with red spray head and white cap, internal ascending tube with slide valve, protection seal



Contents:	purified non-perishable water (gamma irradiated), 620 ml, Limited good for use period: 2 years, when seal is unbroken only	
Materials:	Bottle, protection cap: Tubes (inside and outside): Spray head Tipping valve body Tipping valve ring:	Polyethylene Polyethylene Polyethylene Polystyrene Polyoxymethylene
Dimensions:	approx. 90 x 70 x 250 mm	
Weight:	approx. 750 grams	
Marking:	**BARIKOS KS** – EN 15154-4 - **620 ml** – Charge „xx“ (Charge-No.) – „**yyyy/mm**“ (end of “good for use”-time) – **954004** (Article-No.)	
Handling/Storage:	At room temperature! Protect from cold, heat and direct solar radiation!	
Warning advice:	Eye rinsing does not replace an ophthalmological examination	

Stand 10/2019

BartelsRieger Atemschutztechnik GmbH
Richard-Byrd-Straße 23 · 50829 Köln
Telefon +49 221 59777-0 · Fax +49 221 59777-159
www.bartels-rieger.de · E-Mail: mail@bartels-rieger.de

Die Augenspülflasche