Profile™ 3000

Nitrile Gloves



MASTER CODE	SPECIFICATIONS	PACKAGING	PRODUCT CERTIFICATIONS
PNA	240mm (9.5") Nitrile	100 pieces per dispenser box; 10 dispenser boxes per carton; Recyclable Components	EN 421:2010 · ISO 374-5: 2016 VIRUS · ISO 374-1/ Type C · CE 0598

UK e-mail: info@nitritex.com
Tel: +44 (0) 1638 663338
Fax: +44 (0) 1638 668890
Nitritex Ltd, CB8 7YY, United Kingdom

Nitritex
www.nitritex.com

Profile™ 3000

Nitrile Gloves

MASTER CODE


PNA

RE-ORDER CODE	SIZE
* PNA3100-XXS	XX-Small
PNA3101-XS	X-Small
PNA3102-S	Small
PNA3103-M	Medium
PNA3104-L	Large
PNA3105-XL	X-Large
PNA3106-XXL	XX-Large
PNA3107-XXXL	XXX-Large

SINGLE WALL THICKNESS MINIMUM (mm)

CUFF: 0.07
PALM: 0.09
FINGERTIP: 0.11


PRODUCT ILLUSTRATION



White

COLOUR:	White
LENGTH:	240mm (9.5")
MATERIAL:	Nitrile
TYPE:	Ambidextrous
SURFACE:	Textured
OTHER FEATURES:	Powder-free, latex-free, beaded cuff
PROTEIN LEVEL:	Nitrile contains no natural latex proteins
PACKING:	Boxed
AQL:	1.5 freedom from pinholes (1000ml water test)
PHYSICAL PROPERTIES:	Complies with the requirements of European Glove Standards EN420:2003 +A1:2009 and EN374-1:2016, EN374-5:2016 and EN421:2010
QUALITY SYSTEMS:	Manufactured in a facility holding ISO 9001:2015 and (EU)2016/425 Module D certifications, suitable for contact with all foodstuffs. Regulation (EC)1935/2004 and (EC) 2023/2006.
CLASSIFICATION:	Conforms to Category 3 - Complex Design Personal Protective Equipment PPE Regulation (EU)2016/425
SAMPLING & INSPECTION:	In accordance with ISO 2859
STORAGE INSTRUCTIONS:	Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light
BAR CODES:	Each size carries a unique bar code on the inner dispenser box (EAN-13) and the carton (ITF-14) to allow automatic data entry into Stock Control Systems
LOT NUMBER:	This is shown on the inner dispenser box and the carton to provide product traceability
PACKAGING:	100 pieces per dispenser box; 10 dispenser boxes per carton (1000 pieces)
NOTES:	*Please call for minimum order quantity

PACKAGING ILLUSTRATION



Contacts:

Europe: europe@profile-omega.com
Asia: malaysia@profile-omega.com
Website: www.profile-omega.com

Profile™