

ThermoFisher SCIENTIFIC

Product Specifications

Product **406800025**
1,2-Dichloroethane,99+%,pure

General Product Data

Version	4
CAS No	107-06-2
Molecular weight	98.96
Molecular formula	C2 H4 Cl2
Linear formula	ClCH2CH2Cl
Flash point (°C)	13

Product Specifications

Appearance (Form)	Clear liquid
Color scale	=<50 APHA
Infrared spectrum	Conforms
Refractive index	1.4440 to 1.4460 (20°C, 589 nm)
Water	=<0.05 %
GC	>=99.0 %

Product Specifications are subject to amendment and may change over time.

Acros Organics BV

ENA23, zone 1, nr 1350, Janssen Pharmaceuticalaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: http://www.thermofisher.com
1 Reagent Lane, Fair Lawn, NJ 07410,USA
Fax 201-796-1329