Acetone

AR

Material No: 6776-45
Effective from: 2013/12/05

Product Specification

Test	Specification
Assay	>= 99.5 %
Water (KF)	<= 0.3 %
Titrable Acid (meq/g)	<= 0.0005
Residue after Evaporation	<= 0.001 %
Alcohols	<= 0.06 %
Substances insoluble in water	passes test
Substances Reducing KMnO4 (as O)	<= 0.0002 %

For Laboratory, Research or Manufacturing use.

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000