Dimethyl Sulfoxide
ChromAR® HPLC

Material No: 2969-25
Effective from: 2013/06/29

Product Specification

Test	Specification
Assay (GC)	>= 99.5 %
Water (KF)	<= 0.05 %
Titrable Acid (meq/g)	<= 0.0005
Residue after Evaporation	<= 0.001 %
Ultraviolet Absorbance (1.00-cm cell vs water)	
270 nm	<= 0.4
280 nm	<= 0.2
300 nm	<= 0.1
320 nm	<= 0.05
330 nm	<= 0.02
350 nm	<= 0.01

For Laboratory, Research or Manufacturing use.
Suitable for Liquid Chromatography and UV-Spectrophotometry

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000