

ThermoFisher SCIENTIFIC

Product Specifications

Product **152880250**
Eosin Y,pure,high purity, biological stain

General Product Data

Version	2
CAS No	17372-87-1
Molecular weight	691.85
Molecular formula	C20 H6 Br4 Na2 O5
Linear formula	

Product Specifications

Appearance (Color)	Red to brown
Appearance (Form)	Crystalline powder
Infrared spectrum	Conforms
Dye content	>=90 % (on dry substance)
Loss on drying	=<10 %

Product Specifications are subject to amendment and may change over time.

Acros Organics BV

ENA23, zone 1, nr 1350, Janssen Pharmaceuticalaan 3a, B-2440 Geel, Belgium
Tel +32 14/57.52.11 - Fax +32 14/59.34.34 Internet: http://www.thermofisher.com
1 Reagent Lane, Fair Lawn, NJ 07410,USA
Fax 201-796-1329