Cyclohexane
ChromAR HPLC

Material No: 6720-25
Effective from: 2025/12/19

Product Specification

Test	Specification
Assay (GC)	>= 99.5 %
Residue after Evaporation	<= 0.0005 %
Titrable Acid (meq/g)	<= 0.0005
Water (KF)	<= 0.01 %
Ultraviolet Absorbance (1.00-cm cell vs water)	
210 nm	<= 0.6
220 nm	<= 0.3
230 nm	<= 0.08
240 nm	<= 0.04
350 nm	<= 0.004

For Laboratory, Research or Manufacturing use.
Suitable for Liquid Chromatography and UV-Spectrophotometry

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000