Florisil

Baker Analyzed

Activated

Material No: 7061-00
Effective from: 2015/09/15

Product Specification

Test	Specification
Florisil recovery test with stearyl-stearate/mineral oil standard acc. to ISO 9377-2	passes test
Physical Data (not specifications):	
Average Particle Diameter, um (APD)	140 - 250
Bulk Density (g/cc)	0.5

For Laboratory, Research or Manufacturing use.
For Mineral Oil Determination

Barbara Rykowska
Quality Assurance Coordinator

For questions on this Product Specification please contact Technical Service at technical.service@avantorsciences.com
Avantor Performance Materials Poland S.A. Sowinskiego 11 Gliwice 44-101 Poland Tel: +48 322392000