Seit 1927
über 90 Jahre

Since 1927
more than 90 Years

**Amarell**
GmbH & Co. KG
Präzisions-Thermometer / Aräometer

Lindenstraße 43
D-97892 Kreuzwertheim

**Tel.** +49 (0) 9342 9283-0
**Fax:** +49 (0) 9342 39860

**Internet:** www.amarell.de
**E-Mail:** thermometer@amarell.de

# Technical Data Sheet

| | |
|---|---|
| Art. No.: | A301860-FL |
| Type: | Thermometer similar to ASTM 130 C, prismatic, white backed |
| Specification: | Tank |
| Measurement range: | -7+105°C |
| Subdivision: | 0,5°C |
| Length: | 303mm |
| Diameter: | 6-7mm |
| Filling: | Red special liquid |
| Immersion: | Total |
| Certificate: | Without |

September'2021