In order to provide you with the best possible website experience, we are planning system maintenance from Friday, February 11, 8:00 PM to Saturday, February 12, 9:00 AM CST (Saturday, February 12, 3:00 to 16:00 CET). During this time certain site features and functionality (Order Placement, Pricing and Availability and Quality Documents) will not be available. We know this may interrupt your scheduled work and apologize for any inconvenience. If you do experience issues with site features and functionality, we kindly ask that you check back once updates are complete.

Specification Sheet

Product Name	Methyl heptadecanoate analytical standard
Product Number	51633
Product Brand	SIAL
CAS Number	1731-92-6
Molecular Weight	284.48

TEST	SPECIFICATION
APPEARANCE (COLOR)	White to Off White
APPEARANCE (FORM)	Powder or Crystals
REM. ON APPEARANCE	MAY ALSO BE SOLIDIFIED MASS

TEST	SPECIFICATION
PURITY (GC AREA %)	≥ 99.0 %
PROTON NMR SPECTRUM	CONFORMS TO STRUCTURE
EXPIRY PERIOD	36 MONTHS

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.
Reproduction of any materials from the site is strictly forbidden without permission.
Site Use Terms | **Privacy Policy** | **General Terms and Conditions of Sale** | **Copyright Consent**