PRODUCT DATA SHEET

OMNICAN 20 0.5 ML INS U-40 30GX5/16" (9161619S)



Attribute Name	Attribute Value
Brand name	Omnican®
Product name (HC)	Omnican® 20
Product type	Syringe set
Application	Insulin
Construction of syringe	3-piece
Nominal capacity	0.5 ml
Total graduated capacity	0.5 ml
Max. usable capacity	0.50 ml
Scale value	0.0125 ml
Scale value (I.U.)	0.50 IE
Concentration (device)	20 IE
Light-protected	No

B. Braun Melsungen AG
Carl-Braun-Straße 1
34212 Melsungen | Germany
Tel. +49 5661 710
www.bbraun.com

Attribute Name	Attribute Value
Color of light protection	-
with Cannula	integrated
Outer-ø cannula	0.30 mm
Gauge (cannula)	30
Length of cannula (metric)	8.00 mm
Length of cannula (inch)	5/16"
Type of filter	No filter
Position of filter	-
Type of packaging	Blister
Lubricants	Silicone oil
Packaging contains PVC	No
Packaging contains Latex	No
Packaging contains DEHP	No
Product (unpack.) contains PVC	No
Product (unpack.) cont. Latex	No
Product (unpack.) cont. DEHP	No
Product contains PVC	No
Product contains Latex	No
Product contains DEHP	No
Product contains Nickel	Yes
Remark critical attributes	Needle cont. Nickel
Sterility	sterile
Type of Sterilization	EO
Disposable / Reusable	disposable